
Trash Page 1 of 12





COLLOIDAL SILVER

YOUR IMMUNE SYSTEM'S BEST FRIEND

Colloidal Silver, which is a solution in which microscopic, electrically charged particles of Silver are in suspension, is known to kill 650 different kinds of pathogens, such as bacteria, fungus, viruses, streptococcus and many others, by inhibiting the enzyme which allow their oxygen metabolism. And it does so in minutes, safely, with no side effect.


PRODUCE YOUR OWN COLLOIDAL SILVER AT NO COST

Colloidal Silver can cost as much as 30 dollars for a bottle with only a few ounces. You can now make your own for the price of water.

PROTECT AND IMPROVE THE HEALTH OF YOUR FAMILY, FRIENDS, PETS AND PLANTS.

Trash Page 2 of 12





Many conditions have benefited from the use of Colloidal Silver, such as gastritis, malaria, parasitic infections, psoriasis, eye, ear and mouth infections, yeast infections, either with internal or external use, Colloidal Silver has been as well very successful in veterinary use, for instance the dreaded canine parvo virus. Or you can use a diluted solution either to spray on plant leaves or to treat the soil. And you can purify a gallon of water in minutes by adding 1/2 teaspoon of Colloidal Silver.



ALWAYS HAVE A RELIABLE, FRESH SOURCE AT YOUR FINGERTIPS, ANYTIME, ANYWHERE, FOR YEARS.

IT IS QUICK:

IT IS EASY TO OPERATE:

Just plug the cable into the generator's output, connect the alligator clips to the Silver wires, wait a few minutes, it's done!

YOU WILL GET YEARS OF USE:

For a start, enough Pure Silver to last you for the next two years.

Colloidal Silver -

The Rediscovery of a Super Antibiotic?

Colloidal silver appears to be a powerful, natural antibiotic and preventative against infections. Acting as a catalyst, it reportedly disables the enzyme that one-celled bacteria, viruses and fungi need for their oxygen metabolism. They suffocate without corresponding harm occurring to human enzymes or parts of

Trash Page 3 of 12

the human body chemistry. The result is the destruction of disease-causing organisms in the body and in the food.

Early Research

Colloidal silver was in common use until 1938. Many remember their grandparents putting silver dollars in milk to prolong its freshness at room temperature. At the turn of the century, scientists had discovered that the body's most important fluids are colloidal in nature: suspended ultra-fine particles. Blood, for example, carries nutrition and oxygen to the body cells. This led to studies with colloidal silver. Prior to 1938, colloidal silver was used by physicians as a mainstream antibiotic treatment and was considered quite "high-tech." Production methods, however, were costly. The pharmaceutical industry moved in, causing colloidal research to be set aside in favor of fast working and financially lucrative drugs.

The Food and Drug Administration today classifies colloidal silver as a pre-1938 drug. A letter from the FDA dated 9/13/91 states: "These products may continue to be marketed . . . as long as they are advertised and labeled for the same use as in 1938 and as long as they are manufactured in the original manner." Some of the manufacturing methods used before 1938 are still used today. An electro-colloidal process, which is known to be the best method, is used.

Contemporary Studies

While studying regeneration of limbs, spinal cords and organs in the late 1970s, Robert O. Becker, M.D., author of The Body Electric, discovered that silver ions promote bone growth and kill surrounding bacteria. The March 1978 issue of Science Digest, in an article, "Our Mightiest Germ Fighter," reported: "Thanks to eye-opening research, silver is emerging as a wonder of modern medicine. An antibiotic kills perhaps a half-dozen different disease organisms, but silver kills some 650. Resistant strains fail to develop. Moreover, silver is virtually nontoxic." The article ended with a quote by Dr. Harry Margraf, a biochemist and pioneering silver researcher who worked with the late Carl Moyer, M.D., chairman of Washington University's Department of Surgery in the 1970s: "Silver is the best all-around germ fighter we have."



How It Works

The presence of colloidal silver near a virus, fungus, bacterium or any other single celled pathogen disables its oxygen metabolism enzyme, its chemical lung, so to say. Within a few minutes, the pathogen suffocates and dies, and is cleared out of the body by the immune, lymphatic and elimination systems. Unlike pharmaceutical antibiotics, which destroy beneficial enzymes, colloidal silver leaves these tissue-cell enzymes intact, as they are radically different from the enzymes of primitive single-celled life. Thus colloidal silver is absolutely

Trash Page 4 of 12

safe for humans, reptiles, plants and all multi-celled living matter.

Product Quality

The highest grade is produced by the electro-colloidal / non-chemical method where the silver particles and water have been colloided, i.e., dispersed within and bound to each other by an electric current. The super-fine silver particles are suspended indefinitely in demineralized water. The ideal color of colloidal silver is a golden yellow. The container and dropper must be glass, as plastic cannot preserve the silver in liquid suspension for any length of time.

High concentrations of silver do not kill disease germs more effectively than the safe range of 3 to 5 parts per million (ppm.).

Ingesting Colloidal Silver

Taken orally, the silver solution is absorbed from the mouth into the bloodstream, then transported quickly to the body cells. Swishing the solution under the tongue briefly before swallowing may result in faster absorption. In three to four days the silver may accumulate in the tissues sufficiently for benefits to begin. Colloidal silver is eliminated by the kidneys, lymph system and bowel after several weeks. If routinely exposed to dangerous pathogenic germs, some recommend a regular daily intake as a protection. In cases of minor burns, an accumulation of colloidal silver may hasten healing, reducing the possibility of scar tissue and infection. It is believed by many in the natural healing arts that the lives of millions of people who are susceptible to chronic low-grade infections can be enhanced by this preventative health measure.

Chronic or Serious Conditions

1 teaspoon of 5 ppm. colloidal silver equals about 25 micrograms (mcg.) of silver. 1 - 4 teaspoons per day (25 - 100 mcg.) is generally considered to be a "nutritional amount" and is reported to be safe to use for extended periods of time. Amounts higher than this are generally considered "therapeutic amounts" and should only be used periodically.

In cases of illness, natural health practitioners have often recommended taking double or triple the "nutritional amount" for 30 to 45 days, then dropping down to a smaller maintenance dose. Amounts from 1 - 32 ounces per day have reportedly been used in acute conditions.


If your body is extremely ill or toxic, do not be in a hurry to clear up everything at once. If pathogens are killed off too quickly, the body's five eliminatory channels (liver, kidneys, skin, lungs and bowel) may be temporarily overloaded, causing flu-like conditions, headache, extreme fatigue, dizziness, nausea or aching muscles. Ease off on the colloidal silver to a smaller amount and increase your distilled water intake. Regular bowel movements are a must in order to relieve the discomforts of detoxification. Resolve to reduce sugar and saturated fats from the diet, and exercise more. Given the opportunity, the body's natural ability to heal may amaze you.

Topical Uses

Some have used colloidal silver in a nasal spray mister - to reach the sinuses and

Trash Page 5 of 12

nasal passages. Spray bottles have been used for topical use on kitchen and bathroom surfaces, skin, sore throat, eyes, burns, etc. Colloidal silver is painless on cuts, abrasions, in open wounds, in the nostrils for a stuffy nose, and even in a baby's eyes because, unlike some antiseptics, it does not destroy tissue cells. It's excellent as an underarm deodorant, since most underarm odor is caused by bacteria breaking down substances released by the sweat glands!



Some Common Uses of Colloidal Silver

Natural health practitioners have for years recommended taking one tablespoon daily, for four days, to establish a level, then one teaspoon daily for maintenance (proportional to body weight for children). After six weeks, a pause of several weeks has also been recommended by some natural healing arts doctors. Also, colloidal silver can be applied directly to cuts, scrapes, and open sores, or on a bandage for warts. It can be applied on eczema, itches, acne or bug bites. To purify water, add one tablespoon per gallon, shake well and wait six minutes. Mixed this way, it's tasteless. It is not an allopathic poison.

Veterinary and Garden
Use
Colloidal silver has
worked just as well on
pets of all kinds. Used in
proportion to body
weight, it should bring

Trash Page 6 of 12



the same results. In the garden, field or greenhouse, add enough to the water or soil - and the plants will do the rest.



Tolerance To Disease Organisms

We have all heard of the "super-germs" that are resistant to most modern antibiotics. Some believe that single-celled germs cannot mutate into silver-resistant forms, as happens with conventional antibiotics. Therefore no tolerance to colloidal silver would develop through mutation. Also, colloidal silver has not been demonstrated to interact or interfere with other medicines being taken. Inside the body, silver apparently does not form toxic compounds or react with anything other than a germ's oxygen-metabolizing enzyme. Colloidal silver may truly be a safe, natural remedy for many of mankind's ills.



Additionally, there has never been a drug interaction reported between colloidal silver and any other medication. It's difficult to overdose - unless large amounts are ingested. Colloidal silver has been reported by users to be both a remedy and a prevention for numerous infections, colds, flus, and fermentations due to various bacteria, viruses or fungi, even the non-apparent low-grade, general body infections many people have. Living organisms are in the colloidal chemical state, not the crystalline state. Substances already in that form may be more readily assimilated by the body. Colloidal silver is the most useable form of a reputedly effective germ fighter.

A colloidal suspension is ultra-fine particles of one substance, suspended by an electric charge in another substance. Homogenized milk and aerosol sprays are colloidal suspensions. Colloidal silver is pure, metallic silver (not a chemical compound) of particles 15 atoms or fewer, each with a positive electric charge, and attached to a molecule of simple protein. This new particle floats in pure water. The electric charge is stronger than gravity so the silver particles don't sink.

Colloidal Silver in Advance of Illness?

When the possibility of germ exposure is higher, colloidal silver can be taken orally each day or applied topically when there is a skin problem. It's like having

Trash Page 7 of 12

a second defense system. The silver acts only as a catalyst and is stabilized. It is non-toxic, except to one-cell plants and animals, and is non-addicting. It also apparently kills parasites because they have a one cell egg stage in their reproductive cycle.

Older folks reportedly feel younger because their body energies are used for other uses than constantly fighting disease. Digestion has also been reportedly better. Medical research has shown that silver promotes more rapid healing, with less scar tissue, even in the case of severe burns. Successes have been reported in cases that previously have been given up by established doctors. Colloidal silver is tasteless and won't sting even a baby's eyes, and won't upset your stomach.

More Than 650 Diseases Helped?

Colloidal silver has been reported to kill 650 micro-organisms, many of which are associated with human diseases. This does not automatically mean that taking colloidal silver will "cure" diseases "caused" by these germs. Colloidal silver only kills micro-organisms when they are in contact with it for a sufficient period of time. The human body is a complex system which may prevent high enough concentrations of colloidal silver from reaching the "affected area".



The basic guideline that has been recommended for using colloidal silver is that it usually "works" if you can get a high enough concentration to the "affected area". Some will want to experiment with "higher amounts" (such as 8 or more ounces at a time) to find out what it takes to accomplish this. Do not use colloidal silver if you are allergic to contact with silver metals, or if you notice any digestive upset after use.

The following is a partial list of the more than 650 diseases that colloidal silver has been reputed to be successful against: acne, AIDS (Reference 8), allergies, appendicitis, arthritis, athlete's foot, bladder inflammation, blood parasites, blood poisoning, boils, burns, cancer (References 2, 4, 7), candida, cholera, colitis, conjunctivitis, cystitis, dermatitis, diabetes (Reference 1), dysentery, eczema, fibrositis, gastritis, gonorrhea, hay fever, herpes, impetigo, indigestion, keratitis, leprosy, leukemia, lupus, lymphangitis, Lyme disease, malaria, meningitis, neurasthenia, parasitic infections: viral, fungal and bacterial pneumonia, pleurisy, prostate, pruritus ani, psoriasis, purulent opthalmia, rhinitis, rheumatism, ringworm, scarlet fever, septic conditions of the eyes, ears,

Trash Page 8 of 12

mouth, and throat, seborrhea, septicemia, shingles, skin cancer, staphylococcus and streptococcus infections, stomach flu, syphilis, thyroid, tuberculosis, tonsillitis, toxemia, trachoma, all forms of virus, warts, whooping cough, yeast infection, stomach ulcer, canine parovirus and other veterinary uses, and fungal and viral attacks on plants. Simply spray diluted silver on the leaves and add to the soil.

FUNCTION CONTROLS

MODE

Continuous Power (cont) - allows you to have the unit operate over a longer period of time and allows you to manually time the process (best for GOLDEN YELLOW solutions).

Timer - allows you the choice of three pre-timed operations.

TIMER

Set the MODE switch to timer then select the required time 5min,15min, 30min.

The unit will automatically switch off after the desired time. (best for SILVER/GREY solution)

CURRENT

This function is used mainly when making the GOLDEN YELLOW solution. It regulates the amount of current being sent to the silver wires.

Different currents will make different sized clusters of silver ions (colloidal silver).

We have found that a Colloidal Silver solution made up of the different sizes is more likely to be effective against a larger range of pathogens.

To achieve the best results make up the three different solutions and mix them together at the end of the process.

INDICATOR LIGHTS

Power light appears when Unit is plugged into the Mains Supply.

OutPut Status light confirms current is being supplied to the Silver Wires (The Unit is generating)

TO MAKE A GOLDEN YELLOW SOLUTION

Trash Page 9 of 12

Always use distilled water

DAY ONE

Select CONT (CONTINUOUS) POWER in the MODE settings

SOLUTION ONE

For two to five litres

Place silver wires into water container, taking care not to let them touch each other.

DO NOT ALLOW THE ALLIGATOR CLIPS TO BE IMMERSED IN THE DISTILLED WATER, as this will cause them to rust.

Select CURRENT at 50MA

To begin the process, press the RED START button.

Allow the unit to generate for approximately 90 minutes. The OUTPUT STATUS light confirms that power is being fed to the silver wires.

To interrupt the process, or at the conclusion of the desired time, push the RESET button.

Finally, FILTER the water using a fine filter to remove any solids and place in a dark cupboard overnight.

The solution should turn a pale yellow colour by the next day.

SOLUTION TWO

For two to five litres

Select CURRENT at 100MA Repeat process as above

SOLUTION THREE

For two to five litres

Select CURRENT at 200MA Repeat process as above

DAY TWO

Filter all 3 solutions again using a fine filter to remove any solid particles.

Mix all three solutions together.

Do not attempt to make all three solutions, at once, in the same vessel.

The solution is now ready for use.

We have found that a Colloidal Silver solution made up with the three solutions is more likely to be effective against a larger range of pathogens.

Trash Page 10 of 12

TO MAKE A SILVER GREY SOLUTION

Always use distilled water

Select TIMER on the MODE switch

SELECT DESIRED TIME SETTING ON THE TIMER SWITCH

To make 250mls select 5mins To make 500mls select 15mins To make one litre select 30mins

CURRENT

For 5 or 15 min solutions use the 100MA setting

For 30mins solution use the 200mA setting

These solutions are ready for immediate use.

SUGGESTED APPLICATIONS

The most efficient use of either solution is as an aid to strengthening the immune system.

These solutions can be used as both a preventative and in acute situations

GOLDEN YELLOW SOLUTION

A few drops directly under the tongue as an Immune stimulant weekly. In acute symptoms, a few drops may be used 2 to 3 times daily until symptoms subside.

Place 5mls into 250mls of water to use as a disinfectant for cleaning wounds, abrasions and burns (can be applied undiluted if required).

Add 50mls to one litre of water and use as a disinfectant wash for floors etc or use directly from a Spray bottle onto surface.

SILVER GREY SOLUTION

250mls once a week for immune system stimulant.

250mls 2 to 3 times daily until symptoms subside.

Can be used undiluted as a disinfectant for cleaning wounds, abrasions and burns.

Golden Yellow Solution or Silver Grey Solution?


250mls of the 5 minute silver/grey solution is equal in dosage to 5mls of the golden yellow solution.

However, it must be noted that by making a golden yellow solution using the 3 separate current outputs and finally mixing them together, this should be more

Trash Page 11 of 12

effective against a greater range of pathogens.

Your choice of solution may be dependent on the time you have to make it.



Become an Altered States Affiliate

DISCLAIMER

Altered States products are sold for learning, self-improvement and simple relaxation. No statement contained in this catalogue, and no information provided by any Altered States employee, should be construed as a claim or representation that these products are intended for use in the diagnosis, cure, mitigation, treatment or prevention of disease or any other medical condition. The information contained in this catalogue is deemed to be based on reliable and authoritative report. However, certain persons considered experts may disagree with one or more of the statements contained here. Altered States assumes no liability or risk involved in the use of the products described here. We make no warranty, expressed or implied, other than that the material conforms to applicable standard specifications.

The publisher does not accept any responsibility for the accuracy of the information or the consequences arising from the application, use, or misuse of any of the information contained herein, including any injury and/or damage to any person or property as a matter of product liability, negligence, or otherwise. No warranty, expressed or implied, is made in regard to the contents of this material. No claims or endorsements are made for any drugs or compounds currently marketed or in investigative use. This material is not intended as a guide to self-medication. The reader is advised to discuss the information provided here with a

Trash Page 12 of 12

doctor, pharmacist, nurse, or other authorized healthcare practitioner and to check product information (including package inserts) regarding dosage, precautions, warnings, interactions, and contraindications before administering any drug, herb, radionics tool, or supplement discussed herein.